

## Recommended Shelf Life

### ROMIL-BiO™ BioPure Solvents & Reagents

| Code  | Description   | Shelf Life | Opened   |
|-------|---|------------|----------|
| H051  | Acetonitrile 10 BiO   | 2 years    |          |
| H053  | Acetonitrile 30 BiO   | 2 years    |          |
| H055  | Acetonitrile 50 BiO   | 2 years    |          |
| D511  | Acetonitrile anhydrous BiO  | 2 years    | 3 months |
| D325  | Acetonitrile BiO with DTT   | 2 years    |          |
| D315  | Acetonitrile 10% BiO aqueous solution                             | 1 year     | 6 months |
| D330  | Acetonitrile 20% BiO aqueous solution                             | 1 year     | 6 months |
| D854  | Activator reagent DCI 3% BiO 4,5-dicyanoimidazole/acetonitrile    | 2 years    | 3 months |
| D820  | Activator reagent ETT 3.25% BiO 5-ethylthiotetrazole/acetonitrile | 2 years    | 3 months |
| D513  | Activator reagent TZ 1% BiO tetrazole/acetonitrile                | 2 years    | 3 months |
| D542  | Activator reagent TZ 2.5% BiO tetrazole/acetonitrile              | 2 years    | 3 months |
| D627  | Activator reagent TZ 3% BiO tetrazole/acetonitrile                | 2 years    | 3 months |
| D735  | Activator reagent TZ 3.5% BiO tetrazole/acetonitrile              | 2 years    | 3 months |
| H060  | Ammonia solution BiO  | 2 years    |          |
| D373  | Buffer concentrate BiO acetate/hexanesulphonate                   | 2 years    | 6 months |
| D519  | Capping reagent A BiO acetic anhydride/lutidine/thf               | 2 years    | 3 months |
| D556  | Capping reagent A BiO acetic anhydride/pyridine/thf               | 2 years    | 3 months |
| D610  | Capping reagent A BiO acetic anhydride/thf                        | 2 years    | 3 months |
| D633  | Capping reagent B BiO N-methylimidazole/pyridine/thf              | 2 years    | 3 months |
| D525  | Capping reagent B BiO N-methylimidazole/thf                       | 2 years    | 3 months |
| D318  | 1-Chlorobutane BiO  | 2 years    |          |
| D1682 | Detritylation reagent 3% BiO dca/dichloromethane                  | 2 years    | 3 months |
| D1623 | Detritylation reagent 3% BiO dca/1,2-dichloroethane               | 2 years    | 3 months |
| D1629 | Detritylation reagent 9% BiO dca/1,2-dichloroethane               | 2 years    | 3 months |
| D1754 | Detritylation reagent 1% BiO tca/1,2-dichloroethane               | 2 years    | 3 months |
| D535  | Detritylation reagent 3% BiO tca/dichloromethane                  | 2 years    | 3 months |
| Q535  | Detritylation reagent 3% BiO tca/dichloromethane                  | 2 years    | 3 months |
| H186  | 1,2-Dichloroethane BiO  | 2 years    | 3 months |
| H203  | Dichloromethane BiO   | 2 years    | 3 months |
| H251  | Dimethylformamide BiO   | 1 year     | 3 months |
| H254  | Dimethylformamide BiO with molecular sieve                        | 2 years    | 3 months |
| H240  | Di-iso-propylethylamine BiO                                       | 2 years    |          |
| D102  | Dithiothreitol BiO  | 3 years    |          |
| D312  | Ethyl Acetate BiO   | 2 years    | 3 months |
| H351  | Formamide BiO   | 2 years    | 3 months |
| D306  | n-Heptane BiO   | 2 years    |          |
| H440  | 3-Methylbutan-1-ol BiO  | 2 years    |          |
| D354  | N-Methylpiperidine solution BiO butanol/water/propan-2-ol         | 1 year     | 3 months |
| D332  | N-Methylpiperidine solution BiO methanol/water                    | 1 year     | 3 months |
| H567  | N-Methyl-2-pyrrolidone BiO  | 2 years    | 3 months |
| D643  | Oxidiser reagent BiO iodine in water/lutidine/thf                 | 2 years    | 3 months |
| D532  | Oxidiser reagent BiO iodine 0.02M in water/pyridine/thf           | 2 years    | 3 months |
| D701  | Oxidiser reagent BiO iodine 0.1M in water/pyridine/thf            | 2 years    | 3 months |
| H620  | Phenol BiO glass distilled  | 2 years    |          |
| D342  | Phenylisothiocyanate 5% BiO n-heptane solution                    | 1 year     | 1 month  |
| H613  | Piperidine BiO  | 2 years    | 3 months |
| D501  | Piperidine 20% BiO dmf solution                                   | 1 year     | 3 months |
| H649  | Pyridine BiO  | 1 year     | 3 months |
| H651  | Pyridine BiO with molecular sieve                                 | 2 years    | 3 months |
| D370  | Tetrahydrofuran 3.5% BiO aqueous solution                         | 1 year     | 3 months |
| D368  | Tetrahydrofuran 5% BiO aqueous solution                           | 1 year     | 3 months |
| D369  | Tetrahydrofuran 10% BiO aqueous solution                          | 1 year     | 3 months |
| D349  | Trifluoroacetic Acid BiO  | 2 years    | 3 months |
| D502  | Trifluoroacetic Acid 3% BiO aqueous solution                      | 1 year     | 3 months |
| D362  | Trifluoroacetic Acid 25% BiO aqueous solution                     | 1 year     | 3 months |
| D336  | Trimethylamine 12.5% BiO aqueous solution                         | 1 year     | 3 months |
| D503  | Wash solution BiO acetic acid/water                               | 2 years    |          |

*Recommended Shelf Life*

*ROMIL-BiO™ BioPure Solvents & Reagents*

| Code | Description                                 | Shelf Life | Opened |
|------|---|------------|--------|
| D506 | Waste neutraliser BiO ethanolamine/methanol | 3 years    |        |

For more information see [www.romil.com/stab.htm](http://www.romil.com/stab.htm)

Published: 31 Oct 2007