

Recommended Shelf Life

ROMIL PUROM™ High Purity Process Solvents

Code	Description	Shelf Life	Opened
P5014	Acetic Acid PUROM	2 years	3 months
P5031	Acetone PUROM	2 years	3 months
P5048	Acetonitrile PUROM	2 years	
P5046	Acetonitrile PUROM for preparative HPLC	2 years	
P5083	Butan-1-ol PUROM	2 years	
P5140	Chloroform PUROM stabilised with amylene	2 years	
P5202	Dichloromethane PUROM	2 years	3 months
P5220	Diethyl Ether PUROM stabilised with BHT	2 years	3 months
P5253	Dimethylformamide PUROM	2 years	3 months
P5297	1,4-Dioxan PUROM	2 years	3 months
P5314	Ethanol absolute PUROM	2 years	
P5346	Ethyl Acetate PUROM	2 years	
P5367	n-Heptane 95% PUROM	2 years	
P5366	n-Heptane 99% PUROM	2 years	
P5390	Hexane fraction PUROM	2 years	
P5388	iso-Hexane 95% PUROM	2 years	
P5389	n-Hexane 95% PUROM	2 years	
P5410	Methanol PUROM	2 years	
P5447	Methyl tert-Butyl Ether PUROM	2 years	3 months
P5446	Methyl iso-Butyl Ketone PUROM	2 years	
P5569	n-Octane 95% PUROM	2 years	
P5571	n-Pentane 95% PUROM	2 years	
P5601	Petroleum Ether 40-60° C PUROM	2 years	
P5602	Petroleum Ether 60-80° C PUROM	2 years	
P5624	Propan-1-ol PUROM	2 years	
P5625	Propan-2-ol PUROM	2 years	
P5702	Tetrachloroethylene PUROM	2 years	
P5718	Tetrahydrofuran PUROM	2 years	3 months
P5719	Tetrahydrofuran PUROM stabilised with BHT	2 years	3 months
P5771	Toluene PUROM	2 years	
P5901	2,2,4-Trimethylpentane PUROM	2 years	
P5950	Water PUROM	2 years	
P5982	Xylene mixed isomers PUROM	2 years	

For more information see www.romil.com/stab.htm

Published: 14 Jul 2008