

Recommended Shelf Life

ROMIL-SpS™ Super Purity Solvent Mixes

Code	Description	Shelf Life	Opened
R3502	GC Extraction mix CE1 SpS cyclohexane/ethyl acetate 50:50	1 year	
R4035	HPLC Eluant mix AA1 SpS acetonitrile with acetic acid 0.1%	1 year	
R7651	HPLC Eluant mix AB1 SpS acetonitrile/buffer 35:65	1 year	
R7662	HPLC Eluant mix AB2 SpS acetonitrile/buffer 40:60	1 year	
R7663	HPLC Eluant mix AB3 SpS acetonitrile/buffer 60:40	1 year	
R2579	HPLC Eluant mix AF1 SpS acetonitrile with formic acid 0.1%	1 year	
R4703	HPLC Eluant mix AP1 SpS acetonitrile/propan-2-ol 88:12	1 year	
R0902	HPLC Eluant mix AT1 SpS acetonitrile with trifluoroacetic acid 0.1%	1 year	
R4647	HPLC Eluant mix AT2 SpS acetonitrile with trifluoroacetic acid 0.02%	1 year	
R3033	HPLC Eluant mix AT3 SpS acetonitrile with trifluoroacetic acid 0.05%	1 year	
R7623	HPLC Eluant mix AW1 SpS acetonitrile/water 50:50	1 year	
R4912	HPLC Eluant mix AW2 SpS acetonitrile/water 60:40	1 year	
R6017	HPLC Eluant mix AW4 SpS acetonitrile/water 70:30	1 year	
R8239	HPLC Eluant mix AW5 SpS acetonitrile/water 30:70	1 year	
R7850	HPLC Eluant mix AWA1 SpS acetonitrile/water 95:5 with ammonia 0.1%	1 year	
R6254	HPLC Eluant mix AWF2 SpS acetonitrile/water 5:95 with formic acid 0.1%	1 year	
R5339	HPLC Eluant mix AWF3 SpS acetonitrile/water 95:5 with formic acid 0.1%	1 year	
R5168	HPLC Eluant mix AWP1 SpS acetonitrile/water 37:63 with sodium perchlorate 2.4%	1 year	
R5179	HPLC Eluant mix AWP2 SpS acetonitrile/water 58:42 with sodium perchlorate 2.4%	1 year	
R7620	HPLC Eluant mix MB1 SpS methanol/buffer 25:75 pH 3.2	1 year	
R7684	HPLC Eluant mix MB2 SpS methanol/buffer 25:75 pH 3.75	1 year	
R5215	HPLC Eluant mix MB3 SpS methanol/buffer 20:80 with tetrahydrofuran 1.5%	1 year	
R4076	HPLC Eluant mix MF1 SpS methanol / buffer	1 year	
R4618	HPLC Eluant mix MF2 SpS methanol/water 90:10 / buffer	1 year	
R8927	HPLC Eluant mix MF3 SpS methanol/water 10:90 / buffer	1 year	
R3891	HPLC Eluant mix MF6 SpS methanol with formic acid 0.1%	1 year	
R2013	HPLC Eluant mix MT1 SpS methanol with trifluoroacetic acid 0.1%	1 year	
R5226	HPLC Eluant mix MTB4 SpS methanol/tetrahydrofuran/buffer 69:14:17	1 year	
R4145	HPLC Eluant mix MW1 SpS methanol/water 30:70 with acetone 0.5%	1 year	
R4132	HPLC Eluant mix MW2 SpS methanol/water 30:70	1 year	
R7636	HPLC Eluant mix MW3 SpS methanol/water 50:50	1 year	
R7185	HPLC Eluant mix MWP1 SpS methanol/water 65:35 with phosphoric acid 0.6%	1 year	
R6712	HPLC Eluant mix TPA1 SpS trimethylpentane/propan-2-ol 98:2 with acetic acid 0.1%	1 year	
R6126	HPLC Eluant mix WA2 SpS water with ammonia 0.1%	1 year	3 months
R2803	HPLC Eluant mix WAC1 SpS water with ammonium hydrogen carbonate 10mM pH 9	1 year	3 months
R3680	HPLC Eluant mix WF1 SpS water with formic acid 0.1%	1 year	3 months
R1803	HPLC Eluant mix WT1 SpS water with trifluoroacetic acid 0.1%	1 year	3 months
R2416	HPLC Eluant mix WT2 SpS water with trifluoroacetic acid 0.05%	1 year	3 months
R3527	HPLC Eluant mix WT3 SpS water with trifluoroacetic acid 0.02%	1 year	3 months
R5237	HPLC Wash solution MW4 SpS methanol/water 50:50	1 year	
F3333	MCTS Calibration mix A3 SpS tri-component reference material MCTS elemental analyser	1 year	
R1524	Tetrahydrofuran/Ethanol mix SpS 90:10 by volume	1 year	

Recommended Shelf Life

ROMIL-SpS™ Super Purity Solvent Mixes

Code	Description	Shelf Life	Opened
R3126	Trifluoroacetic Acid/Water 10%/3% SpS in dichloromethane	1 year	
R1121	Trifluoroacetic Acid/Water 15%/3% SpS in dichloromethane	1 year	
R1387	Trifluoroacetic Acid/Water 20%/3% SpS in dichloromethane	1 year	
R2005	Trifluoroacetic Acid/Water 25%/3% SpS in dichloromethane	1 year	
R1303	Trifluoroacetic Acid 1M SpS aqueous solution	1 year	

For more information see www.romil.com/stab.htm

Published: 17 Aug 2008