

Recommended Shelf Life

ROMIL-SpR™ Ion-Pair & Buffer Reagents

Code	Description	Shelf Life	Opened
HR079	Ammonium Acetate SpR	3 years	
IP231	Benzenesulphonic Acid sodium salt SpR	3 years	
IP204	Butane-1-sulphonic Acid sodium salt SpR	3 years	
IP230	(DL)-Camphor-10-sulphonic Acid sodium salt SpR	3 years	
IP232	D-Camphor-10-sulphonic Acid sodium salt SpR	3 years	
IP210	Decane-1-sulphonic Acid sodium salt SpR	3 years	
IP352	Di-n-butylamine SpR	2 years	
IP524	Di-n-butylamine Phosphate SpR 0.5M solution	2 years	
IP226	Dioctyl Sulphosuccinate sodium salt SpR	3 years	
IP212	Dodecane-1-sulphonic Acid sodium salt SpR	3 years	
IP412	Dodecyltrimethylammonium Bromide SpR	3 years	
IP222	Ethane-1,2-disulphonic Acid di-sodium salt SpR	3 years	
IP235	Heptafluorobutyric Acid SpR	2 years	
IP207	Heptane-1-sulphonic Acid sodium salt SpR	3 years	
IP216	Hexadecane-1-sulphonic Acid sodium salt SpR	3 years	
IP416	Hexadecyltrimethylammonium Bromide SpR	3 years	
IP206	Hexane-1-sulphonic Acid sodium salt SpR	3 years	
IP8C	Ion-Pair Kit SpR for CATIONIC separations	3 years	
IP9A	Ion-Pair Kit SpR for ANIONIC separations	3 years	
IP228	Lithium Dodecyl Sulphate SpR	3 years	
IP201	Methanesulphonic Acid sodium salt SpR	3 years	
IP224	Naphthalene-2-sulphonic Acid sodium salt SpR	3 years	
IP209	Nonane-1-sulphonic Acid sodium salt SpR	3 years	
IP208	Octane-1-sulphonic Acid sodium salt SpR	3 years	
IP205	Pentane-1-sulphonic Acid sodium salt SpR	3 years	
HA622	Perchloric Acid 60% SpR	2 years	
HR706	Potassium Bromide SpR for infra-red spectroscopy	2 years	
IP227	Potassium Ethyl Sulphate SpR	3 years	
HR471	di-Potassium Hydrogen Phosphate 3H2O SpR	3 years	
IP203	Propane-1-sulphonic Acid sodium salt SpR	3 years	
IP223	Propane-2-sulphonic Acid sodium salt SpR	3 years	
IP229	Sodium Dodecyl Sulphate SpR	3 years	
IP404	Tetrabutylammonium Bromide SpR	3 years	
IP504	Tetrabutylammonium Hydroxide SpR 0.5M solution	2 years	
IP518	Tetrabutylammonium Hydroxide SpR 40% solution	2 years	
IP454	Tetrabutylammonium Phosphate SpR	3 years	
IP554	Tetrabutylammonium Phosphate SpR 0.5M solution	2 years	
IP558	Tetrabutylammonium Phosphate SpR 1M solution	2 years	
IP414	Tetradecyltrimethylammonium Bromide SpR	3 years	
IP402	Tetraethylammonium Bromide SpR	3 years	
IP452	Tetraethylammonium Hydrogen Sulphate SpR	3 years	
IP401	Tetramethylammonium Chloride SpR	3 years	
IP451	Tetramethylammonium Hydrogen Sulphate SpR	3 years	
IP408	Tetraoctylammonium Bromide SpR	3 years	
IP403	Tetrapropylammonium Bromide SpR	3 years	
IP354	Tri-n-butylamine SpR	2 years	
IP356	Tri-n-hexylamine SpR	2 years	
IP358	Tri-n-octylamine SpR	2 years	
IP355	Tri-n-pentylamine SpR	2 years	

For more information see www.romil.com/stab.htm

Published: 31 Oct 2007