

Recommended Shelf Life

ROMIL-UpS™ Ultra Purity Solvents

Code	Description	Shelf Life	Opened
H033	Acetone UpS ultra residue	2 years	3 months
H050	Acetonitrile UpS ultra gradient	2 years	
H054	Acetonitrile UpS ultra residue	2 years	
H204	Dichloromethane UpS ultra gradient	1 year	3 months
H205	Dichloromethane UpS ultra residue	1 year	3 months
H317	Ethanol absolute UpS ultra gradient	2 years	
H347	Ethyl Acetate UpS ultra gradient	2 years	
H364	n-Heptane 99% UpS ultra residue	2 years	
H391	n-Hexane 95% UpS ultra gradient	2 years	
H394	n-Hexane 99% UpS ultra residue	2 years	
H412	Methanol UpS purge and trap	2 years	
H411	Methanol UpS ultra gradient	2 years	
H415	Methanol UpS ultra residue	2 years	
H573	n-Pentane 99% UpS ultra residue	2 years	
H626	Propan-2-ol UpS ultra gradient	2 years	
H720	Tetrahydrofuran UpS ultra gradient	2 years	3 months
H772	Toluene UpS ultra residue	2 years	
H903	2,2,4-Trimethylpentane UpS ultra residue	2 years	
H949	Water UpS ultra gradient	2 years	1 month

For more information see www.romil.com/stab.htm

Published: 14 Jul 2008